


LAB MANUAL: EXPERIMENT 2

Aim: To verify Lambert-Beer's law using a given solution of potassium dichromate at the wavelength of its maximum absorption (λ_{\max}) and consequent determination of the unknown concentration of a solution of potassium dichromate.

Theory:

The Beer-Lambert law states that the absorbance of a solution is directly proportional to the concentration of the absorbing species in the solution and the path length. Thus, for a fixed path length (cuvette length), UV/Vis spectroscopy can be used to determine the concentration of the absorber in a solution. The absorbance changes with concentration. Thus, a higher concentration of the colored solution absorbs more light (and transmits less) than a solution of lower concentration.




$$\log (I_0 / I_t) = A = \epsilon c l$$

According to Beer-Lambert law,


$\log(I_0/I_t)=A=\epsilon cl$ where I_0 and I_t are the incident and transmitted intensities,

A = absorbance and ϵ is a constant i.e. absorptivity (also called the extinction coefficient).

If the concentration is measured in molL^{-1} , the absorptivity is called molar absorptivity. $A= \epsilon cl$. At constant length $A \propto c$



Working principle of spectrophotometer



Spectrophotometer



Cuvette

Requirements:

spectrophotometer, cuvette, six test tubes, Measuring cylinder, 10 mL pipette, 0.001M $K_2Cr_2O_7$ solution, distilled water, test tube rack, and tissues (preferably lint-free).

Procedure & observation table:

Step 1: To record the absorbance of $K_2Cr_2O_7$ solution at different wavelengths to determine the light wavelength for its maximum absorption (λ_{max}):

- (a) Prepare 200 mL of 0.001M $K_2Cr_2O_7$ (**Molecular weight 294.18 gm/mol**) solution in distilled water.
- (b) Label five clean, dry, test tubes 1–5.
- (c) Use a 10 mL pipette to prepare five standard solutions according to Table 1.
- (d) Thoroughly mix each solution.
- (e) Calibrate the spectrophotometer with respect to the blank solution i.e. distilled water.
- (f) Fill the first one of the prepared solutions (1-5) up to a certain level in the cuvette of the spectrophotometer.
- (g) Record the absorbance of the respective solution at different wavelengths as mentioned in Table 2.
- (h) Plot the absorbance data in the graph paper with respect to the wavelength and calculate the light wavelength for its maximum absorption (λ_{max}) in $K_2Cr_2O_7$.

Table 1:

Test-tube	0.001M $K_2Cr_2O_7$ (mL)	Distilled water (mL)	Concentration (M)
1	1	9	
2	2	8	
3	3	7	
4	4	6	
5	5	5	

Table 2: The solution of the **No.** the test tube was chosen for the determination of the light wavelength for its maximum absorption (λ_{\max}).

Entry	Wavelength (λ in nm)	Absorbance
1		
2		
3		
4		
5		
6		

Step 2: To record the absorbance of different concentrations of solutions at the specified λ_{\max} :

- Set the operating wavelength of the spectrophotometer in the range of absorption maxima of aqueous $K_2Cr_2O_7$ solution (λ_{\max}).
- Calibrate the spectrophotometer with respect to water as the blank.
- Fill each of the solutions up to a certain level in the cuvette of the spectrophotometer.
- Record the absorbance of the respective solutions as stated in Table 3.
- Plot the absorbance data in the graph paper with respect to the concentration which should be a straight line passing through the origin.

Table 3:

Entry	Test-tube	Absorbance
1	1	
2	2	
3	3	
4	4	
5	5	

Step 3: Determination of the unknown concentration of a given potassium dichromate solution:

- (a) Fill the solution up to a certain level in the cuvette of the spectrophotometer.
- (b) Record the absorbance of the given solution of unknown concentration.
- (c) Plot the absorbance data in the same graph obtained above (ideally it should be on the same straight line obtained from the plot of step 1)
- (d) Draw a perpendicular line from the absorbance point to the concentration axis.
- (e) Note down the corresponding unknown concentration.

Conclusion:

1) The light wavelength for its maximum absorption (λ_{\max}) is found to be **nm.**

2) The concentration of the unknown solution was found to be **M**

Precautions:

- (a) Always mix the standard solutions properly.
- (b) Wipe the outside of the cuvette every time with a lint-free tissue.
- (b) Handle cuvettes only by the top edge of the ribbed sides.
- (c) Dislodge any bubbles by gently tapping the cuvette on a hard surface.
- (d) Always position the cuvette so the light passes through the clear sides.
- (e) Always set the light source of the instrument in the absorption maxima range of the given solution.